



OPTICAL SOFTWARE  
USER MANUAL (DRAFT)

REFERENCE	QU 5659
REVISION AND ISSUE:	Rev A
DATE:	9th April 2018
CLIENT:	Imperial College
AUTHOR:	Mike Parkin
BUSINESS DEVELOPMENT MANAGER:	Chris Garratt



<b>Document Summary</b>	<b>3</b>
<b>Software Overview</b>	<b>4</b>
Starting the Application	4
User Interface - Overview	7
User Interface - Scope	9
User Interface - DAQ	9
User Interface - Monochromator	10
User Interface - Shutter	10
User Interface - Bottom Pane	11
Advanced Configuration	12
Executing a Test - Kinetic Mode	14
Executing a Test - Spectrum Mode	15
Stopping a Test	16
File I/O - Data Files	16
File I/O - Error Log	17
<b>Code Overview</b>	<b>18</b>
LabVIEW Project File	18
Top Level VI Overview	19
Navigating the Architecture	20
Creating a New Module	21
Adding a New Module	30
Removing/Disabling a Module	40
Testing an Isolated Module	41
Further Reading	42
Appendix - File I/O Handling	43

## **Document Summary**

This document aims to give an overview on how to use the Optical Software Application and a description of the main screens contained within in it. It also provides a high level description of the LabVIEW architecture that has been used to write it and how modifications can be made to it as additional features and requirements arise in the future.


## Software Overview

### Starting the Application


In order to start the application, the LabVIEW project file must first be opened. This is what links the application to all of its dependencies.

To do this, navigate to where the LabVIEW project file (Optical Software.lvproj) is located and double click to open it.


**NOTE:** This application has been developed in LabVIEW 2017 so must be opened with this version of LabVIEW. It is not possible to open it with a previous version of LabVIEW.



This will open the LabVIEW Project Explorer window as seen below. Double click 'Top Level.vi' to open the window from where the application can be started.




This will open the window seen below. From here the user can set whether to run with real or simulated hardware. Press the Grey arrow to start the application and see the UI.



## User Interface - Overview

When the application is started the user will be presented with a UI as seen below. This is where various test parameters can be entered and where tests can be started and stopped. During a test, data will be visible on the graphs and the indicators on the left hand side will be updated.



The key thing to notice is that anything with a **white** background is a control where information *can* be entered by the user. Anything with a **grey** background is an indicator and *cannot* be updated by the user.

### Control


Data can be entered by the user...

### Indicator

Data cannot be entered by the user...

## User Interface - Left Pane

The pane on the left hand side of the UI contains a mixture of both controls and indicators where the user can set certain test parameters and also observe the progress of ongoing tests.



The screenshot shows a vertical list of controls in a light gray panel. At the top is a 'Folder' label above a text box containing 'C:\Austin Alpha\Optical Software Code Simplified\Output' and a folder icon. Below is a 'Sample Name' label above a text box with '08-05-18'. Then a 'Path to Files' label above a text box with 'C:\Austin Alpha\Optical Software Code Simplified\Output\08-05-18'. This is followed by a 'Mode' section with radio buttons for 'Spectrum' and 'Kinetic' (selected). Below that is a 'B Test?' checkbox which is checked. Then a 'Number of Averages' label above a spinner box set to '12'. Next is an 'Amplification Factor' label above a text box with '400'. Then a 'Calculation Method' label above a dropdown menu showing 'Transmission'. Below that is a '100% Reflectance Standard' label above a text box with '1'. The bottom section is titled 'Current Test' and shows 'A-Test' in a box. Below this are four more text boxes: 'A-Test Counter' (0), 'B-Test Counter' (0), 'Background Level (V)' (0), and 'Background Level (Amps)' (0).

**Folder** - The directory specified by the user where any data files will be saved.

**Sample Name** - The sub-folder within the Folder location where any data files will be saved.

**Path to Files** - The location of the data files for that test (Folder + Sample Name).

**Mode** - Allows the operator to select whether to run tests in Spectrum or Kinetic mode.

**B Test?** - Choose whether to run a B Test in addition to an A Test.

**Number of Averages** - The number of averages/acquisitions per test.

**Amplification Factor** - The value defined by the user that links the Background Level (V) to Background Level (Amps).

**Calculation Method** - The type of calculation that is applied at the end of each test.

100% Reflectance Standard - Value used in the Kubelka-Munk calculation.

**Current Test** - Informs the operator what test is currently running.

**A-Test Counter** - Informs the operator how many averages have currently been acquired for an A-Test.


**B-Test Counter** - Informs the operator how many averages have currently been acquired for a B-Test.

**Background Level (V)** - The background voltage calculated from DAQ data at the end of an A-Test.

**Background Level (Amps)** - The background level in Amps (Background Level (V) x Amplification Factor).

The User Interface displays a selection of tabs showing high level configuration parameters and controls for the devices connected to the system. The ability to modify advanced configurations is shown later in the document.

## User Interface - Scope




The screenshot shows the 'Scope' tab selected in a four-tab interface (Scope, DAQ, Monochromator, Shutter). Below the tabs, there are two input fields with up/down arrows:

- Scope Vertical Range (V) with a value of 0.
- Scope Vertical Offset (0.0 divs) with a value of 0.4.

**Scope Vertical Range (V)** - The vertical range setting sent down to the scope.

**Scope Vertical Offset (0.0 divs)** - The vertical offset setting sent down to the scope.

## User Interface - DAQ



The screenshot shows the 'DAQ' tab selected in a four-tab interface (Scope, DAQ, Monochromator, Shutter). Below the tabs, there are six input fields with up/down arrows:

- Rate (Hz) with a value of 250000.
- dt with a value of 4E-6.
- Average Number of Points with a value of 1.
- Number of Samples with a value of 275000.
- End Time (s) with a value of 1.
- Final Number of Points with a value of 275000.

**Rate (Hz)** - The rate at which the DAQ will sample.



**Average Number of Points** - The number of points that the acquired DAQ data will be averaged by.

**End Time** - The length of time the DAQ will sample for.

**$\Delta t$**  - The time between each data point acquired.

**Number of Samples** - Total number of samples ((Rate x End Time) x 1.1).

**Final Number of Points** - Total number of points after averaging.

## User Interface - Monochromator

**Current Wavelength (nm)** - The current wavelength that the monochromator is/will be showing.

**Desired Wavelength (nm)** - The wavelength, entered by the operator, that the monochromator will be moved to when the **Move Monochromator** button is pressed.

**Move Monochromator** - Moves the monochromator to the wavelength specified

## User Interface - Shutter

**Move Shutter** - Moves the shutter based on settings in the advanced config.

## User Interface - Bottom Pane

At the bottom of the UI there is a series of buttons that the operator uses to Start and Stop tests as well as Exiting the application safely. When a test is started certain buttons on the UI are disabled to prevent the operator mistakenly interfering with the current test.



**Current Config** - Displays the configuration currently being used by the application.

**Start Test** - Starts and runs a test based on the configuration set by the operator

**Stop Test** - Stops the current test that is running.

**Advanced Config** - Takes the operator through the Advanced Configuration screen which specified low level settings of the devices connected to the system - See *Advanced Configuration* section.


**Exit** - Clears the buffers and Exits the application safely.

## Advanced Configuration

When the application is started it will *always* load the default configuration file. All settings and configurations are saved to a config file. The default configuration file is saved to the same location as the project file and cannot be saved or overwritten from the application. However, it can be used as a basis for operators to create and save their own configuration file.

Selecting the **Advanced Config** button enables the user to make changes to the configuration settings currently being displayed, save it as a new file or load an existing configuration file.

The operator can set a configuration that suits and save it using the **Save Config** button. If the operator tries to overwrite the default configuration file, a pop up will warn them to save it as a new one with a different name.



The application will then guide the operator through the process of saving the new file - which will then be used by the application and can then be used in the future.

To load an existing configuration, select the **Load Config** button. The operator will select the configuration file and populate the application with the settings contained in that file. This must be done every time the application is started.



The Advanced Configuration screen allows the operator to modify all settings that aren't shown on the main UI. These settings include things such as the COM Port associated with a device or the type of Scope and Monochromator being used in the system - things that are unlikely to change frequently.

Whenever a setting is changed the operator must **Save Config** in order for the new settings to take effect.

# Imperial College London

FileIO
Scope
DAQ
Monochromator
Colour Wheel
Shutter

VISA Resource USB0:

Model DPO 3012

Timebase 0.0004

Timeout (ms) 30000

**Trigger**

Source Channel 2

Threshold (V) 0.25

Mode Laser (TAS)

**PIA mode**

LED line 1

Test LEDs ☐

LED on (s) 1000

LED off (s) 1000

**TAS**

Channel Channel 1

Vertical Coupling DC

Probe Attenuation (10.0) 1

Vertical Range (V) 0

Vertical Offset (0.0 divs) 0.4

Bandwidth FULL

Impedance 50

Multiplication factor 1

**TPC**

Channel None

Vertical Coupling DC

Probe Attenuation (10.0) 1

Vertical Range (V) 0

Vertical Offset (0.0 divs) 0.1

Bandwidth FULL

Impedance 1M

Multiplication factor 0

**Scope Indicators**

**Model**

Acquisition Type Hi Res

Record length (points) 10,000

Averages on scope 1

Average number of points 1

Final number of points 10000

Pretrigger points 1000

Current Config: C:\Project Work\icl-optical-software-rewrite\Optical Software Code Simplified\Config.ini

Load Config

Save Config

Exit

Selecting the **Exit** button will take the operator back to the main User Interface.

## Executing a Test - Kinetic Mode

If the test has been started in Kinetic mode, the operator will be asked to confirm the current wavelength displayed by the monochromator and the desired wavelength of the monochromator. I.e, the wavelength will this test be run at.

When confirmed, the monochromator will move to the requested wavelength and the colour wheel will automatically move to the appropriate setting. The DAQ and Scope will then start acquiring based on the trigger mechanism specified in the configuration.

During the test, the UI will display the number of tests completed (for either A-Test and B-Test).

Current Test	A-Test
A-Test Counter	7
B-Test Counter	0
Background Level (V)	0
Background Level (Amps)	0


At the end of an A-Test, the data acquired will be processed and data files generated at the location specified by the operator. The filename of these data files will be prefixed with the wavelength that they were testing at. By this point the Background Level (V) and Background Level (Amps) will have been calculated and updated on the UI.

The application will then automatically move the shutter and start a B-Test. The B-Test, if uninterrupted, will always run the same number of averages as the A-Test; meaning the counters should always match at the end (the screenshot below is during a B-Test).

Current Test	B-Test
A-Test Counter	10
B-Test Counter	2
Background Level (V)	0.062808
Background Level (Amps)	25.1232

At the end of a B-Test, the data acquired will be processed in the same way as the A-Test. However, this time there will also be A-B calculation stage and an AB Stitching stage. These will be written to a file and saved to the same location.

The shutter will then return to its original position and a pop up will appear informing the operator that the acquisition has completed.




## Executing a Test - Spectrum Mode

Spectrum Mode is similar in the fact that it runs A-Tests and B-Tests across a number of averages as specified by the operator.



However, the main difference is that it runs multiple series of A and B Tests consecutively, each time at a different wavelengths. Eg, A-Test at 600nm, B-Test at 600nm, Data processed and files created, A Test at 700nm, B-Test at 700nm, Data files processed etc...

When the Mode on the UI is specified as Spectrum and a test is started, the operator will be asked to enter the current wavelength of the monochromator. This is followed by a request to define the desired wavelengths. The operator must set the **Start Wavelength**, **End Wavelength** and **Wavelength Interval** and press OK. A series of wavelengths is created and shuffled into a random order that will be used to determine what wavelength the monochromator and colour wheel moves to at the start of each test.



The UI will give an indication as to what wavelength is currently being tested and how far through that test it is. There is also the also the option to check the data file directory to see which files have been processed.

The process, once started, requires little operator interaction as stages of the test will happen automatically until there are no more wavelengths to be tested. The Acquisition Complete message will be displayed when all tests have completed.

## Stopping a Test

Stopping the test can be done by selecting the **Stop Test** button on the UI. This will often be the only button enabled whilst a test is running and can be pressed at any point regardless of what test is being run.

If the current test is an A-Test, pressing the **Stop Test** button will stop the A-Test and move onto the B-Test. The B-Test will then run the same number of averages that the A-Test reached when it was stopped. The post-processing calculations will then be executed as normal.

The Stop Test button is not limited to which test it can stop and can stop both tests if needed.




## File I/O - Data Files

Throughout the testing procedure the raw and final data files of both the DAQ and Scope are generated and written to the location specified by the operator. There is also a parameters file that lists the key parameters set during this test - this is different to the application configuration file.

During a test there will be some temporary files located here. This is for performance reasons and means data that will be used at a later stage in the test is saved to a file rather than the application holding it in memory. In the event of a PC failure, this data will not be lost and could provide explanation as to any fault.

An overview of the File I/O from a coding perspective can be found in the Appendix of this document.



**File I/O - Error Log**

The application has a dedicated Error Module that will write any errors to an error log. The location of this is defined in the Advanced Configuration screen. It is not possible to view this log through the application itself but the file can be opened by notepad.



## Code Overview

### LabVIEW Project File


All the files used in this application are now contained within a LabVIEW project file (.lvproj). Utilising a LabVIEW project file provides a way of organising and managing an application as it grows in size and complexity.

As an application grows, the number of modules and SubVIs will also increase. Having these VI's and dependencies collated in an easy to read format means any developers or users involved in this project know where to look should modifications need to be made.

This environment can be used to open the Top Level VI and any module contained within it. In the case of this application, each module is responsible for a certain area of functionality.

Each module has its own LabVIEW Library (.lvlib). The benefit of using a LabVIEW Library is that all files associated with that module are contained within that library, from the APIs (used for sending messages to this module from other modules) to SubVI's (blocks of code that have been wrapped to contain certain functionality or purpose).

In the image to the right, it can be seen that there is a 'UI' library. This contains 'Main.VI'. This is the VI that contains the SubVI's and is the receiver of any messages sent to it using the API's.




In this application the DAQ module sends information to the UI module. Therefore, 'Set DAQ data API.vi' is used to send this data. If we wanted to send a confirmation from the UI module back to the DAQ module, we would need to use an API from the DAQ library.

In addition to this being a more intuitive layout for the developer it also provides benefits when it comes to debugging. Should an issue arise relating to a certain device in the system, it can be easily identified where, in the application, this issue is likely to be originating from. For example, if there is an issue with DAQ data being read it is highly probable the problem lies in the DAQ Module. From here, various debug techniques can be used to deduce what the problem is.

More information and features of the Project Explorer can be found in the following white paper on the National Instruments website: <http://www.ni.com/white-paper/7197/en/>.

## Top Level VI Overview

The application is run from the Top Level VI. This must be opened from Project Explorer and can be seen below. The block diagram of the Top Level VI is seen in below.



This provides an indication of how the architecture has been split into modules. Each module runs in parallel and provides a different purpose. Existing modules include: Main Controller Module, User Interface Module, Shutter Module, Monochromator Module, Colour Wheel Module, DAQ Module, Scope Module, File IO Module and Post Processing Module.


Breaking it down like this is better for performance and is good practice. Rather than one loop being responsible for all functionality, there are multiple loops sharing the work.

It should be made clear that the Top Level VI **is not** the User Interface. The User Interface is a module in itself and is the only module that is configured to show its Front Panel on execution. All other modules (and the Top Level VI) will run in the background and will not show their Front Panels.

## Navigating the Architecture

The layout of each module will look similar to an extent. That is based on the fact that they all originate from the same 'Example Module'. This has been left in the project file and the 'Creating a Module' section will explain how it is possible to use this to create new modules.

The block diagram of the Example Module is shown below.



The module uses 'User Events' to interact and transfer messages with other modules. The event structure has dynamically registered for an event and is constantly waiting for a message to be sent to it. As mentioned previously, this message will have been sent using one of the API's from within the Example Module library file. When a message is received by the module it unpacks it and reads the command and data. This command is wired directly to a case structure and will perform whatever task is within that case.

This is a highly expandable method of module development in that many different commands/API's can be created and the module will be able to handle them; as long as a case has been created within that module.

A module could have dozens of different commands that it needs to handle. It is down to the developer as to what this command will do when the module receives it.


More information on Event-Driven Programming in LabVIEW can be found at the following link: <http://www.ni.com/white-paper/3331/en/> and additional information on Dynamically Registering for Events can be found here: [http://zone.ni.com/reference/en-XX/help/371361K-01/lvhowto/dynamic\\_register\\_event/](http://zone.ni.com/reference/en-XX/help/371361K-01/lvhowto/dynamic_register_event/)

## Creating a New Module


When looking to create a new module it is best to make a duplication of the 'Example Module' library. As it is a LabVIEW Library, all of its contents, layout and structure will be replicated and there will be no name conflicts. This is another benefit of using libraries.

The steps of creating a new module using Project Explorer are listed below:


1. Create a copy of 'Example Module.lvlib' by right clicking on it and selecting Save As.




## 2. Select a new Name and Select Continue...




## 3. Select a directory location and save it.




4. Check that the folder and its contents have been saved correctly.



5. Switch to Project Explorer to view the new library in the project file.




6. For organisation purposes, move this library to the 'Modules' folder by clicking it and dragging.






7. It is sensible to change the appearance of all VI's within this library. Ie, change the color theme of the icons so sub VI's from this module are easily distinguishable when reviewing the code. In order to do this, right click on the library and select Properties.
8. Go to the General Settings category and select Edit Icon...







- Using the toolbox palette on the right hand side, change the color of the banner (at the top of the icon) and use the Icon Text tab to change the text that is displayed. Below is an example. It can be whatever theme suits the developer. The important thing is that it is different from the existing modules. When happy with the theme, select OK.




- A preview of how the VI's in this library will appear is now shown. Select OK apply the changes.




11. A warning will be shown informing that this Icon will be applied to all the VI's in this library and overwrite the existing Icons. Select Yes to confirm this is understood.




12. As these changes affect all VI's in the library, they will need to be saved. When prompted to do this select Save - All.




13. From the Project Explorer, navigate to Main.vi within the new library and double click to open it.



14. In the top right of the window it can be seen that the VI now displays the Icon create in the previous steps.




15. Pressing Ctrl + E will display the Block Diagram of this VI. This is the shell that will receive messages sent from other modules and execute functionality that is contained within it.




## Adding a New Module

16. Register Module for Events. It is now important that we Register this Module so that it can receive messages sent to it. This was explained in more detail earlier on in the document. Navigate to the Events library in the Project Explorer and Open Init Ref Kill.vi. As the name suggests this vi is responsible for initialising, referencing and killing the references used for passing the event based messages around the application.




17. Press Ctrl + E to view the Block Diagram of the VI.



18. Right Click on the block diagram and select Select a VI from the menu. Navigate to Modules/New Module Folder/Controls/Command Data [cluster].ctl and select OK.




19. This will place one of the controls previous created dynamically when the new library was created. This control is a cluster specific to the library and contains the data that will be passed around within every message that is sent. The cluster consists of an enum and variant. This is generic enough to be adaptable to any message and data that needs to be sent.




20. Right click on the Block Diagram and using the search function search for Create User Event. When found, drag this onto the Block Diagram. Alternatively, copy and paste one of the Create User Event functions that is already on the Block Diagram. Wire the New Module control to the User Event Data Type input on the Create User Event function.

21. Right click on the User Event Output of the Create User Event function and create a control.




22. A control will be placed on the Block Diagram.







23. Double click this to find it on the Front Panel.




24. This now needs to be added to the Type Definition that contains all of the existing Event References. On the Front Panel, select the New Module reference and press Ctrl + X to add it to the clipboard. Right click on the Type Definition containing all the existing Event References and select Open Type Def. Press Ctrl + V to paste the New Module reference onto the Front Panel of this Type Definition.




25. Drag the new reference into the cluster containing all the existing Event References. Save this and close the window.



26. Go back to the Block Diagram of the Init Ref Kill.vi. The bundle function will now have an additional input that requires a wire.




27. Link the Event Reference output from the Create New User function to this input.




28. Save this VI and close the window.


29. Go back to Main.vi that is within the New Module library.
30. Switch to the User Event case of the Event Structure and right click on the drop down menu of it. Select Edit Events Handled by This Case...




31. Change the User Event being received by this module by selecting <New Module>: User Event from the Event Sources list. Select OK to apply this.



32. Save Main.vi and close the window.
33. From the Project Explorer, open Top Level.vi. This displays all the Modules currently used by the application.
34. In order to place the New Module in the Top Level vi, right click on the Block Diagram and select Select a VI from the menu. Navigate to Modules/New Module Folder/Main.vi and Select Ok.




35. As with the existing Modules, wire the error cluster output of the Start Events vi to the error cluster input of the New Module Main.vi. Wire the error cluster output of the New Module Main.vi to the Merge Errors node. Expand the Merge Errors if necessary.




## Removing/Disabling a Module

In order to remove a module it can simply be deleted from the Block Diagram of the code. It won't break the code but, from that point on, any API calls made to that module will do nothing and get no response.

A more subtle, and easily reversible, method would be to disable that section of the code. Using quick drop (Ctrl + Spacebar) and searching for Diagram Disable Structure will provide a structure that can be dragged around any part of the code that the developer wants to be able to view but not be included during runtime.



Right-clicking on the Diagram Disable Structure and selecting 'Enable this Subdiagram' will bring the code back into play and be included during runtime.



Right-clicking on the Diagram Disable Structure and selecting 'Disable this Subdiagram' will revert it back.

More information on the Diagram Disable Structure can be found here:

[http://zone.ni.com/reference/en-XX/help/371361P-01/glang/diagram\\_disable\\_structure/](http://zone.ni.com/reference/en-XX/help/371361P-01/glang/diagram_disable_structure/)

## Testing an Isolated Module

The modularity of the architecture means that each module can be run on its own. I.e, if one aspect of the system needs to be tested it can be done so by select the module from the project explorer and executing the code. This should be treated as a debugging tool to debug a certain but of the application rather than a permanent method of controlling one of the devices.

A key thing to remember when doing this is that the config is no longer being sent from the File IO module. Therefore, the developer will need to input this in themselves.



## **Further Reading**

Debugging Techniques in LabVIEW: <http://www.ni.com/tutorial/14130/en/>

Event-Driven Programming in LabVIEW: <http://www.ni.com/white-paper/3331/en/>

Dynamically Registering for Events:

[http://zone.ni.com/reference/en-XX/help/371361K-01/lvhowto/dynamic\\_register\\_event/](http://zone.ni.com/reference/en-XX/help/371361K-01/lvhowto/dynamic_register_event/)

Diagram Disable Structure:

[http://zone.ni.com/reference/en-XX/help/371361P-01/glang/diagram\\_disable\\_structure/](http://zone.ni.com/reference/en-XX/help/371361P-01/glang/diagram_disable_structure/)


## Appendix - File I/O Handling

### A TEST

DAQ




SCOPE



### B TEST

SCOPE




Temporary Data Files

Permanent Data Files


Process Steps

## POST PROCESSING

DAQ



Scope



A-B

